



sutherlandmedical

# Introducing CollarCare Light™





sutherlandmedical

# CollarCare Light™

CollarCare Light™ arm slings are designed to provide support to the clavicle area or to provide full immobilisation of the shoulder.

The CollarCare Light™ arm sling is particularly suitable for CVA patients requiring long term limb support and immobilisation.

The “O” buckle design allows for the CollarCare Light™ arm sling to be adjusted to suit a large variety of applications.

The CollarCare Light™ arm sling also includes a circumferential waist strap allowing the arm to be fixed to the body for complete shoulder immobilisation.

CollarCare Light™ arm slings can be easily fitted by patients without assistance.

CollarCare Light™ arm slings are available in 4 sizes for application on a variety of arm sizes.



## CollarCare Light™

Code	Application	Sling	Size
SM040	Shoulder Immobiliser	16 x 32cm	Small
SM050	Shoulder Immobiliser	19 x 39cm	Medium
SM060	Shoulder Immobiliser	21 x 44cm	Large
SM070	Shoulder Immobiliser	24 x 50cm	Extra Large

### Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria 3166

Phone: 1300 664 027 Fax: 1300 664 028 Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

Or visit us today: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)