



sutherlandmedical

# Introducing CollarForm™ Arm Sling:





sutherlandmedical

# CollarForm™

CollarForm™ is a strong, stretch resistant sling made of a 5cm wide strip of dense foam covered with a pure cotton stockinette.

CollarForm™ is ideally designed to immobilise the clavicle.

CollarForm™ is 100% Australian Made.

Code	Width	Length	Packaging
SM000	5cm	10m	Economy Roll
SM001	5cm	2 x 6m	Dispenser Pack

## Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria

Phone: 1300 664 027 Fax: 1300 664 028 Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

Or visit us today: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)