

# SMV™

## STERILISATION PACKAGING

SMV™ is pleased to offer comprehensive range of premium grade, single use sterilisation packaging for use in public and private hospitals, primary healthcare, veterinary, and other sterilisation applications.

SMV™ sterilisation packaging features premium 70 and 60gsm Arjo Wiggins paper and high quality puncture resistant films with integrated sterilisation indicators on the film face and self seal flap.

SMV™ visual sterilisation packaging is TGA registered and manufactured to EN 11607 standard. The product is accepted and qualified for sale in the USA, Germany and Japan.



## SMV™ Visual Packaging Product Range

- Self Sealing Sterilisation Pouches
- Heat Sealing Sterilisation Pouches (Flat and Gusseted)
- Heat Sealing Sterilisation Reels (Flat and Gusseted)



sutherlandmedical

**Free Call:** 1300 664 027 **Free Fax:** 1300 664 028

**Email:** [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

**Website:** [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)

# SMV™ Visual Packaging Product Range

## Self Seal Products

Code	Description	Qty	Size (mm)
SMV-100	Self Seal Pouch - Flat (60gsm)	Box of 200	57mm x 100mm
SMV-102	Self Seal Pouch - Flat (60gsm)	Box of 200	90mm x 135mm
SMV-103	Self Seal Pouch - Flat (60gsm)	Box of 200	70mm x 230mm
SMV-104	Self Seal Pouch - Flat (60gsm)	Box of 200	90mm x 230mm
SMV-105	Self Seal Pouch - Flat (60gsm)	Box of 200	135mm x 255mm
SMV-106	Self Seal Pouch - Flat (60gsm)	Box of 200	190mm x 330mm
SMV-107	Self Seal Pouch - Flat (60gsm)	Box of 200	250mm x 370mm
SMV-108	Self Seal Pouch - Flat (60gsm)	Box of 200	300mm x 460mm

## Heat Seal Pouches

Code	Description	Qty	Size (mm)
SMV-200	Heat Seal Pouch - Flat (70gsm)	Box of 200	57mm x 130mm
SMV-201	Heat Seal Pouch - Flat (70gsm)	Box of 200	70mm x 250mm
SMV-202	Heat Seal Pouch - Flat (70gsm)	Box of 200	90mm x 250mm
SMV-203	Heat Seal Pouch - Flat (70gsm)	Box of 200	135mm x 275mm
SMV-204	Heat Seal Pouch - Flat (70gsm)	Box of 200	190mm x 350mm
SMV-205	Heat Seal Pouch - Flat (70gsm)	Box of 200	250mm x 370mm
SMV-206	Heat Seal Pouch - Flat (70gsm)	Box of 200	50mm x 200mm
SMV-210	Heat Seal Pouch - Gussetted (70gsm)	Box of 100	100mm x 50mm x 300mm
SMV-211	Heat Seal Pouch - Gussetted (70gsm)	Box of 100	250mm x 60mm x 450mm

## Heat Seal Reels

Code	Description	Qty	Size (mm)
SMV-400	Flat Roll - Heat Seal (70gsm)	Each	55mm x 200m
SMV-401	Flat Roll - Heat Seal (70gsm)	Each	75mm x 200m
SMV-402	Flat Roll - Heat Seal (70gsm)	Each	100mm x 200m
SMV-403	Flat Roll - Heat Seal (70gsm)	Each	150mm x 200m
SMV-404	Gussetted Roll - Heat Seal (70gsm)	Each	100mm x 50mm x 100m
SMV-405	Gussetted Roll - Heat Seal (70gsm)	Each	150mm x 50mm x 100m
SMV-410	Flat Roll - Heat Seal (70gsm)	Each	200mm x 200m
SMV-411	Flat Roll - Heat Seal (70gsm)	Each	300mm x 200m



sutherlandmedical

**Free Call:** 1300 664 027 **Free Fax:** 1300 664 028

**Email:** [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

**Website:** [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)