

Tubular Form Measuring Guide

Limb Circumference

15 - 25cm

26 - 35cm

36 - 44cm

45 - 50cm

51 - 60cm

LOW COMPRESSION

post plaster removal, swelling, oedema
Superficial varices, scalds & burns

Size D

Size E

Size F

Size G

Size H

MODERATE COMPRESSION

post plaster removal, swelling, oedema
varicose conditions, leg ulcer prevention
post burn keloid scarring

Size C

Size D

Size E

Size F

Size G

Note:

All compression is based on a
minimum double layer of Tubular Form

Circumference:

When applying to a limb, measure to
fit the narrowest part.

When applying to the body, measure
the affected area.



sutherlandmedical



Tel: 1300 664 027

E: orders@sutherlandmedical.com.au