



Introducing **THUASNE** orthoses

Premium support, bracing and immobilisation



LOMBAX® HIGH High dorso lumbar support belt

TH0843H40

- 40 cm high
- additional straps with a front-back posture correcting effect, for accurate adjustment of tightness, according to different positions
- aerated Poromax® structure at the back
- two-part front fastening system with hand loops, making the belt easy to put on



LOMBAX® DORSO Posture correcting dorsal lumbar support belt

TH0845

- aluminium structure that conforms to the curve of the spine and can be adjusted to the patient's morphology
- option to adjust the height of the orthosis
- additional straps to adjust compression
- additional strap for sternal support



LOMBATECH® Adjustable reinforced lumbar support belt

TH0851H26

- 26 cm high
- patented, removable, adjustable strap to 32 cm high
- self-adjusting back for optimum anatomical fit to shape
- reinforced abdominal plate for greater support
- easy to put on with hand loop and centering tabs
- patented perspiration-wicking Combirex® fabric

ORTHO CARE



LIGAFLEX® POST-OP OPEN Hinged knee splint [open]

TH2384-OPEN

- 40 to 59 cm high
- adjustable lateral reinforcements
- manual adjustment of flexion from 0° to 120°
- manual adjustment of extension from 0° to 75°
- height adjustable for high-performance support
- single size



LIGAFLEX® POST-OP CLOSED Hinged knee splint [closed]

TH2384-CLOSED

- 50 to 69 cm high
- adjustable lateral reinforcements
- manual adjustment of flexion from 0° to 120°
- manual adjustment of extension from 0° to 75°
- height adjustable for high-performance support
- single size



LOMBAMUM® Posture correcting maternity belt

TH0805

- 21 cm high
- adjustable compression and additional straps for adapted support
- two abdominal support bands so the belt can be worn full term
- abdominal section lined with bamboo fibre fabric for exceptional softness
- two hand loops
- for intermittent and long term use, variable exertion

LOMBAX® HIGH

Code	Size	Measurement
TH0843H40-1	1	52 - 62 cm
TH0843H40-2	2	63 - 74 cm
TH0843H40-3	3	75 - 89 cm
TH0843H40-4	4	90 - 105 cm
TH0843H40-5	5	106 - 120 cm
TH0843H40-6	6	121 - 140 cm

LOMBAX® DORSO

Code	Size	Measurement
TH0845-1	1	52 - 62 cm
TH0845-2	2	63 - 74 cm
TH0845-3	3	75 - 89 cm
TH0845-4	4	90 - 105 cm
TH0845-5	5	106 - 120 cm
TH0845-6	6	121 - 140 cm

LOMBATECH®

Code	Size	Measurement
TH0851H26-1	1	60 - 75 cm
TH0851H26-2	2	76 - 94 cm
TH0851H26-3	3	95 - 115 cm
TH0851H26-4	4	116 - 140 cm

LIGAFLEX® POST-OP OPEN KNEE

Code	Size	Measurement
TH2384-OPEN	Single	40 - 59 cm

LIGAFLEX® POST-OP CLOSED KNEE

Code	Size	Measurement
TH2384-CLOSED	Single	40 - 59 cm

Sutherland Medical Pty Ltd

Phone: 1300 664 027

Fax: 1300 664 028

E-mail: orders@sutherlandmedical.com.au

Web site: www.sutherlandmedical.com.au

LOMBAMUM®

Code	Size	Measurement
TH0805	One	82 - 125 cm