



sutherlandmedical

Introducing **ComfiFast MultiStretch™**
A clever MultiStretch
alternative to conventional
dressing retention:



RETENTION BANDAGES



sutherlandmedical

Comfifast MultiStretch™

Tubular MultiStretch Retention Bandage

Comfifast MultiStretch™ tubular bandages has been developed using MultiStretch technology; enabling the bandage to stretch in multiple directions - radial, longitudinal and diagonal, designed to provide complete freedom of movement and enhance patient comfort.

Comfifast MultiStretch™ is used for:


- Patients with eczema
- Wound Dressing Retention
- Isolation of I.V. Lines and Tubing
- Application to any area where Skin Protection is required for fragile skin

Features:

- Enhanced patient comfort
- Complete freedom of movement
- Comforms to mis-shapen limbs
- Quick and easy to apply
- Fully Washable and cost effective
- Suitable for use under compression bandaging or hosiery garments
- Available in 5 easy to identify colour coded sizes

Code	Size	Application	Roll	Measurement
FMS15	Red	Small Limbs	10m	3.5 cm*
FMS25	Green	Medium Limbs	10m	5 cm*
FMS35	Blue	Large Limbs	10m	7.5 cm*
FMS45	Yellow	X Large Limbs, Head	10m	10.75 cm*
FMS55	Beige	Adult Trunk	10m	17.5 cm*

* Unstretched Width



Fibre content: 92% Viscose, 5% Elastane, 3% Polyamide

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria

Ph: 1300 664 027 Fax: 1300 664 028

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au