


Introducing **MediMaxGlobal**TM Quality disposable ECG electrodes:



Patient
Monitoring

Foam Back ECGs
in three sizes

Resting TAB ECG
Radiolucent ECG



MediMaxGlobal™

Quality disposable ECG electrodes



MediMaxGlobal electrodes are designed for use in emergency, ICU/CCU, surgery and general diagnostic monitoring applications

- Easy to handle for fast application and removal
- Quality Ag/AgCl system providing stable, reliable ECG signals
- Advanced formulation hydrogel for low impedance and improved signal
- Minimal clean up required
- Excellent adhesion whilst hypo-allergenic
- Latex free
- TGA registered

Code	Size	Configuration	Type	Package	Box Qty
MXF55-50	5.5 cm	Round	PE Foam Back	50 per pouch	1,000
MXF43-50	4.3 x 4.5 cm	Oval	PE Foam Back	50 per pouch	1,000
MXF43R-50	4.3 x 4.5 cm	Oval	Radiolucent Foam	50 per pouch	1,000
MXF3238-50	3.2 x 3.8 cm	Oval	PE Foam Back	50 per pouch	1,000
MXT2334	2.3 x 3.4 cm	Card of 10	Resting TAB	100 per pouch	500

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria

Phone: 1300 664 027 Fax: 1300 664 028 Email: orders@sutherlandmedical.com.au

Or visit us today: www.sutherlandmedical.com.au