



sutherlandmedical

Introducing **Tubular-FormSSB™**

The ideal shaped tubular
compression bandage:



COMPRESSION BANDAGES

Tubular-FormSSB™

Tubular Shaped Support Bandage

Tubular-FormSSB™ is ideally designed for the treatment of orthopaedic and vascular conditions..

Tubular-FormSSB™ is used for:

- Acute Venous Insufficiency
- Varicose Ulcers
- Joint Effusion and Soft Tissue Injuries
- Sclerotherapy
- Control of Oedema/Swelling

Features:

- **100% Australian Made**
- Product treated with Unique Low Fray Formula
- Provides Moderate Compression at the ankle and wrist under a single layer
- Easy to use Measuring Guides available
- Unique colour coding system on Full Arm - 1/2 Leg (5 sizes)

Code	Size	Application	Box	Measurement
SM700	Small	Full Arm - 1/2 Leg	10	15 – 19 cm*
SM701	Medium	Full Arm - 1/2 Leg	10	18 – 23 cm*
SM702	Large	Full Arm - 1/2 Leg	10	22 – 27 cm*
SM707	Extra Large	Full Arm - 1/2 Leg	10	26 – 35 cm*
SM708	Extra Extra Large	Full Arm - 1/2 Leg	10	34 – 39 cm*
SM703	Small	Full Leg	10	15 – 19 cm*
SM704	Medium	Full Leg	10	18 – 23 cm*
SM705	Large	Full Leg	10	22 – 27 cm*
SM706	Extra Large	Full Leg	10	26 – 35 cm*



* All measurements are taken at the ankle or wrist.

Fibre Content: Cotton, Elastodiene polyurethane and Polyester

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria

Ph: 1300 664 027 Fax: 1300 664 028

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au