



sutherlandmedical

Introducing **Abdominal Binders:**



2 Panel



3 Panel



4 Panel

SPLINTS AND IMMobilISATION

Abdominal Binders

Abdominal Binders are ideally designed to provide support and compression to the ribs, abdomen or lumbar area in a post-operative or post injury environment.

The Abdominal Binders are manufactured from 7.5cm or 3" panels of heavy grade elastic that are stitched into 2 panel (15cm), 3 panel (22.5cm) and 4 panel (30cm) binders.

The Abdominal Binders have a touch tape closure system with the 15cm suitable for use as a rib binder.

The complete range of Abdominal Binders are 100% Australian Made.

SUPPORT RATING: **FIRM CONTROL**

Application

Measure around the patient's torso in a supine position with a tape measure and select the most appropriate size based on the application guide below. Position the support around the patient with the touch tape closure at the front.

The Sutherland Medical label should be used as a guide to indicate the top of the support. Ensure that the support is firm and comfortable while the patient is standing as well as seated.

Code	Width	Length	Panels	Fitted Waist Measurement
SM830	15cm	Regular	2 panel	75 - 112cm
SM831	15cm	Large	2 panel	105 - 155cm
SM829	22.5cm	Small	3 panel	60 - 90cm
SM832	22.5cm	Regular	3 panel	75 - 112cm
SM833	22.5cm	Large	3 panel	105 - 155cm
SM837	22.5cm	XLarge	3 panel	145 - 170cm
SM834	30cm	Regular	4 panel	75 - 112cm
SM835	30cm	Large	4 panel	105 - 155cm
SM836	30cm	XLarge	4 panel	145 - 170cm

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria

Ph: 1300 664 027 Fax: 1300 664 028

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au