



sutherlandmedical

Introducing CollarCare™ Padded Arm Sling:



SPLINTS AND IMMOBILISATION

CollarCare™

Padded Arm Slings

CollarCare™ padded arm slings are designed to provide full immobilisation of the shoulder and clavicle area.

The CollarCare™ padded arm sling is particularly suitable for CVA patients requiring long term limb support and immobilisation.

The padded fabric provides superior body conformity and weight bearing comfort.

The “D” buckle design allows the CollarCare™ padded arm sling to be adjusted to suit a large variety of applications. The soft padded arm pouch can also be folded to securely fit a variety of limb sizes.

The CollarCare™ padded arm sling also includes a circumferential wrist strap allowing the arm to be fixed to the body for complete shoulder immobilisation. The touchtape velcro closure on the waist strap allows it to be secured anywhere on the fabric.

CollarCare™ padded arm slings can be easily fitted by patients without assistance.

CollarCare™ Padded Arm Slings are made of Polyester, Viscose, Foam, Polycotton, Nylon & Polyamide and is 100% Latex Free and 100% Australian made.

Code	Description	Measurement (Height x Width)	Size
SM991	Clavicle Immobilisation	150mm x 370mm	Child
SM990	Clavicle Immobilisation	210mm x 420mm	Adult

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria
 Ph: 1300 664 027 Fax: 1300 664 028
 Email: orders@sutherlandmedical.com.au
 Web Site: www.sutherlandmedical.com.au