



sutherlandmedical

Introducing **Nitto** Yu-ki Ban™ & Xtrata™ Perme Roll: Surgical tapes & transparent film dressings with STRATEGEL™ technology.



WOUND CARE




Nitto Yu-ki Ban™ & Xtrata™ Perme Roll

Stratagel™ adhesive technology, protecting the skin from both irritation and pain on removal

Above all, Gentle to the skin
Stratum Corneum Protection System

Conventional Product
Skin damage caused by removal of stratum corneum

XTRATA
Very little Stratum corneum is removed



99% of stratum corneum is removed
Conventional Product A

40% of stratum corneum is removed
Xtrata Perme-roll

Features:

XTRATA PERME-ROLL

- Acts as a barrier to external bacteria
- Superior moisture permeability
- Excellent adhesion with low removal of stratum corneum
- Superb flexibility and elasticity to joints and other flexible areas
- Waterproof
- Transparency allows observation of the application site without removal
- Latex Free

YU-KI BAN

- Low irritation tape
- Allows repeated application
- Latex Free



Code	Description	Qty
RL3251	Yu-ki Ban™ Non-woven (water resistant) cloth 12mm x 10m	24 rolls
RL3252	Yu-ki Ban™ Non-woven (water resistant) cloth 25mm x 10m	12 rolls
RL3255	Yu-ki Ban™ Non-woven (water resistant) cloth 50mm x 10m	6 rolls
RL3521	Yu-ki Ban™ Alpha (waterproof) polyolefin film 12mm x 7m	24 rolls
RL3522	Yu-ki Ban™ Alpha (waterproof) polyolefin film 25mm x 7m	12 rolls
RL3525	Yu-ki Ban™ Alpha (waterproof) polyolefin film 50mm x 7m	6 rolls
RLH24R05	Xtrata™ Perme Roll (transparent film dressing) 5cm x 10m	1
RLH24R10	Xtrata™ Perme Roll (transparent film dressing) 10cm x 10m	1

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria
Ph: 1300 664 027 Fax: 1300 664 028
Email: orders@sutherlandmedical.com.au
Web Site: www.sutherlandmedical.com.au