



sutherlandmedical

Introducing **Relevo**: The Super Absorbent Dressing



WOUND CARE

Relevo

The super absorbent dressing for heavily exuding wounds

Relevo dressings absorb 2-3 times more than traditional absorbent dressings.

The design allows the dressing to be applied either side facing the wound.

Relevo absorbs under compression and binds the fluid inside the dressing without leaking.

Features:

- Contributes to fewer changes, time saving and better patient comfort
- Absorbs vertically and can reduce the risk of maceration on healthy skin area
- Contains a combination of cellulose fibres and superabsorbent polymers (SAP) with a polypropylene cover
- Latex free
- Low profile and easy handling

| Code | Description | Absorption capacity | | Qty |
|-------------|-----------------------------|---------------------|------------------------|-----|
| | | (Tap water) | (Saline Solution 0.9%) | |
| RL50011/10 | Relevo Dressing 10cm x 10cm | 120ml | 50ml | 10 |
| RLF50011/50 | Relevo Dressing 10cm x 10cm | 120ml | 50ml | 50 |
| RLF50021 | Relevo Dressing 10cm x 20cm | 300ml | 120ml | 50 |
| RLF50041 | Relevo Dressing 20cm x 20cm | 600ml | 240ml | 50 |



Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria
 Ph: 1300 664 027 Fax: 1300 664 028
 Email: orders@sutherlandmedical.com.au
 Web Site: www.sutherlandmedical.com.au