



sutherlandmedical

# Introducing **Trio** Silvex<sup>®</sup>: silicone convex for soft silicone convexity seal



Feel the difference, experience the benefit of  
our advanced silicone technology.....

OSTOMY CARE



sutherlandmedical

# Silvex<sup>®</sup> silicone convex seal

## Advanced Silicone Technology

Designed with the wearer in mind, **Trio Silvex<sup>®</sup>** offers a soft silicone seal for use when mild convexity is required. The skin-friendly, soft silicone gel allows for greater comfort whilst providing instant adhesion and a perfect custom fit.


**Trio Silvex<sup>®</sup>** will stretch around your stoma filling creases and contours through close conformity, reducing any gaps between the skin and the seal where effluent can escape.

**Trio Silvex<sup>®</sup>** is a soft silicone convex seal that works completely differently to traditional hydrocolloid products. It does not absorb moisture or body waste - for the first time you can have a clean and hygienic seal. **Trio Silvex<sup>®</sup>** convex seals have a natural tendency to return to their original shape when stretched. This means you simply stretch a **Trio Silvex<sup>®</sup>** seal around your stoma and shape it to form a perfect custom fit.

**Trio Silvex<sup>®</sup>** removes in one piece leaving no residue - your pouch change regime can now be quicker, easier and cleaner

### Features:

- Maintains integrity and shape to provide longer lasting seal
- More hygienic than traditional hydrocolloid products - does not absorb faeces or urine
- Allows the skin to breathe
- Reduced risk of leakage - extends wear time of your ostomy product
- Stretch around your stoma and shape to create the perfect custom fit
- **Trio Silvex<sup>®</sup>** conforms to the shape of your stoma
- Promotes comfort due to high level of flexibility
- Instant adhesion - better seal formed
- One piece removal - leaves no residue
- Atraumatic removal - less skin irritation



Code	Description	Qty	SAS Code
TR1020CX	Silvex <sup>®</sup> Silicone Convex Seal, Size 1, 20-30mm	10	80053K
TR1030CX	Silvex <sup>®</sup> Silicone Convex Seal, Size 2, 30-40mm	10	80053K

### Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria  
Ph: 1300 664 027 Fax: 1300 664 028  
Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)  
Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)