




sutherlandmedical

Introducing **Comfast™ Easywrap™**
Garments:
A clever multistretch ready to wear
comfortable garment.



RETENTION BANDAGES



sutherlandmedical

Comfifast™ Easywrap™

Multistretch Retention Garments

Comfifast™ Easywrap™ garments have been developed to help support the treatment of atopic eczema as a skin covering when using steroids, heavy emollients or ointments and for dressing retention eg: burns patients or post surgery.

A full range of garments are available to fit all areas of the body for head to toe dressing. The range has been developed using Multistretch Technology; giving the garment multiple stretch - radial, longitudinal and diagonal, to allow the wearer complete freedom of movement.

Features:

- Acts as a barrier to help prevent the impact of scratching
- Prevents the transfers of emollients, creams and paste to the users' clothes or bedding
- Ready to wear, easy to use and quick to apply. Increases patient compliance with treatment due to simplicity of use, comfort and effectiveness
- Light and breathable fabric allowing the skin to breathe
- Fully washable and can be machine washed up to 30 times for frequent use

PRODUCT	SIZE: LIMB CIRC.	REF
 Adult Vests	Small - Chest: 71 - 91cm (28" - 36") Medium - Chest: 91-117cm (36" - 46") Large - Chest: 117-152cm (46" - 60")	CV20 CV21 CV22
 Adult Leggings	Small - Waist/Hip: 71 - 91cm (28" - 36") Medium - Waist/Hip: 91 - 117cm (36" - 46") Large - Waist/Hip: 117 - 152cm (46" - 60")	CL20 CL21 CL22
 Child Vests	6 - 24 months 2 - 5 years 5 - 8 years 8 - 11 years 11 - 14 years	CV24 CV5 CV8 CV11 CV14
 Child Tights	6 - 24 months	CT24
 Child Leggings	2 - 5 years 5 - 8 years 8 - 11 years 11 - 14 years	CL5 CL8 CL11 CL14
 Child Clavas	6 months - 5 years 5 - 14 years	CC5 CC14
 Child Mittens	Up to 24 months (with thumb hole) 8 years (with thumb hole) 8 - 14 years (without thumb hole)	CM24 CM8 CM14
 Child Socks	Up to 8 years 8 - 14 years	CS8 CS14



Fibre content: 86% Viscose, 11% Nylon, 3% Elastodiene

Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh, Victoria
Ph: 1300 664 027 Fax: 1300 664 028
Email: orders@sutherlandmedical.com.au
Web Site: www.sutherlandmedical.com.au