



sutherlandmedical

Introducing **Abdominal Binders:**



15 cm



22.5 cm



30 cm

SPLINTS AND IMMOBILISATION

Abdominal Binders

Abdominal Binders are ideally designed to provide support and compression to the ribs, abdomen or lumbar area in a post-operative or post injury environment.

Abdominal binders are manufactured from heavy grade elastic and are a fully adjustable one-piece system.

The Abdominal Binders have a hook & loop closure system.

The complete range of Abdominal Binders are 100% Australian Made.

SUPPORT RATING: **FIRM CONTROL**

Application

Measure around the patient's torso in a supine position with a tape measure and select the most appropriate size based on the application guide below. Position the support around the patient with the hook & loop tape closure at the front.

The Sutherland Medical label should be used as a guide to indicate the top of the support. Ensure that the support is firm and comfortable while the patient is standing as well as seated.

Code	Width / Height	Size	Fitted Waist Measurement
SM830*	15cm	Medium	75 - 112cm
SM831*	15cm	Large	110 - 155cm
SM829	22.5cm	Small	60 - 90cm
SM832	22.5cm	Medium	75 - 112cm
SM833	22.5cm	Large	110 - 155cm
SM837	22.5cm	XLarge	145 - 170cm
SM834	30cm	Medium	75 - 112cm
SM835	30cm	Large	110 - 155cm
SM836	30cm	XLarge	145 - 170cm

*The 15cm abdominal binder is also suitable for use as a rib binder.

LATEX

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au