



sutherlandmedical

# Introducing **Avoca** Silver Nitrate Caustic Applicators:



**\*\*\*New \*\*\* Instructions for use**

- If necessary, use **one drop of water**. **Only wet the tip of the applicator for 20 seconds ...no more**

**WOUND CARE**



sutherlandmedical

# Avoca Silver Nitrate Caustic Applicators

## Simple, Effective, Economic

- 75% Silver Nitrate
- Single use 15cm Bendable Applicator
- 100 Applicators in a resealable tube
- Removal of non-genital warts (including verruca's)
- Removal of unwanted tissue
- Cauterisation
- **If necessary, use one drop of water. Only wet the tip of the applicator for 20 seconds ...no more**
- 'how to use' and more technical information available from [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)



Code	Description	Qty
SNA7484	Avoca Silver Nitrate Caustic Applicators	1 tube (100 sticks)

### Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)

ARTG 278284

Version 3 / Created 18/5/2020