



sutherlandmedical

Introducing Cervical Collars:



SM range (SM400S, SM401M, SM402L)



SM Range (SM403, SM404, SM405, SM406, SM407)

SPLINTS AND IMMOBILISATION

Cervical Collars

Soft foam Cervical Collars are designed to limit cervical flexion.

The collars are manufactured from a high density, non latex foam and covered in a natural 100% cotton stockinette to provide a comfortable and contoured fit.

| Code | Size | Application | Height | Length |
|--------|-------------|---------------|--------|--------|
| SM400S | Small | Limit Flexion | 7.5cm | 45cm |
| SM401M | Medium | Limit Flexion | 8.5cm | 50cm |
| SM402L | Large | Limit Flexion | 9cm | 55cm |
| SM403 | Small Long | Limit Flexion | 7.5cm | 53.5cm |
| SM404 | Medium Long | Limit Flexion | 9cm | 53.5cm |
| SM405 | Large Long | Limit Flexion | 10cm | 53.5cm |
| SM406 | Youth | Limit Flexion | 5cm | 52cm |
| SM407 | Child | Limit Flexion | 5cm | 40cm |



Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189
 Ph: 1300 664 027
 Email: orders@sutherlandmedical.com.au
 Web Site: www.sutherlandmedical.com.au