



sutherlandmedical

# Introducing CoFlex<sup>®</sup> TLC:

Two Layer Compression Bandage system designed to meet all your needs.



COMPRESSION BANDAGES

# CoFlex® TLC

## Two Layer Compression

### Non-Medicated

- New itch control first layer. CoFlex® TLC has added a new formula to Layer 1 that now includes aloe to make it more itch-resistant and comfortable on sensitive skin.
- CoFlex® TLC uses an odour control foam 1st layer that absorbs 85% of the odours of standard foam.
- CoFlex® TLC outperforms other 2 layer compression kits from the superior quality of its patented Latex Free cohesive bandage formula and consistent non-constricting compression.

### Medicated

- Layer 1 is a Zinc or Calamine impregnated foam roll that rolls on cleanly and easily as a dry foam roll.
- Layer 2 is the only Latex Free compression bandage to stick directly on a moist dressing allowing the use of a low-profile medicated compression.
- Visual indicators have been incorporated into the cohesive bandage to assist with application in all medicated and non - medicated kits.
- Calamine and Zinc Impregnated Foam rolls are available separately.

#### • For ABI > 0.5

CoFlex TLC Zinc, TLC Calamine, TLC Lite  
25-30mmHg light compression

#### • For ABI > 0.8

CoFlex TLC TN, TLC XL  
35-40mmHg standard compression

Code	Application	Size
A8840UBC-TN	CoFlex® TLC Calamine Lite Kit	L1 - 10cm x 5.5m unstretched, L2 - 10cm x 4.76m unstretched
AUB840C-01	CoFlex® TLC Calamine Medicated Roll	10cm x 5.5m unstretched
A8840UBZ-TN	CoFlex® TLC Zinc Lite Kit	L1 - 10cm x 5.5m unstretched, L2 - 10cm x 4.76m unstretched
AUB840Z-01	CoFlex® TLC Zinc Medicated Roll	10cm x 5.5m unstretched
A7802TLC-TN	CoFlex® TLC Lite Kit	L1 - 10cm x 3.1m unstretched, L2 - 10cm x 3.5m unstretched
A7800TLC-TN	CoFlex® TLC Kit	L1 - 10cm x 3.1m unstretched, L2 - 10cm x 3.5m unstretched
A7800TLC-XL	CoFlex® TLC X/Length Kit	L1 - 10cm x 4.9m unstretched, L2 - 10cm x 4.8m unstretched



#### Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189  
Ph: 1300 664 027  
Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)  
Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)

ARTG 305810 (rolls only)  
ARTG 218947 (Kits only)

Version 4 / Created 18/5/2020