



sutherlandmedical

Introducing **ComfiFAST Multistretch™**:  
A clever multistretch alternative to  
conventional dressing retention.



RETENTION BANDAGES



sutherlandmedical

# Comfifast Multistretch™

## Tubular Multistretch Retention Bandage

**Comfifast Multistretch™** tubular bandages has been developed using Multistretch Technology; enabling the bandage to stretch in multiple tubular directions - radial, longitudinal and diagonal, designed to provide complete freedom of movement and enhance patient comfort.

**Comfifast Multistretch™** is used for:

- Patients with eczema
- Wound Dressing Retention
- Isolation of I.V. Lines and Tubing
- Application to any area where Skin Protection is required for fragile skin

### Features:

- Enhanced patient comfort
- Complete freedom of movement
- Conforms to mis-shapen limbs
- Quick and easy to apply
- Fully Washable and cost effective
- Suitable for use under compression bandaging or hosiery garments
- Available in 5 easy to identify colour coded sizes

Code	Colour	Application	Size	Limb Circumference	Roll
FMS15	Red	Small Limb	3.5cm*	8 – 15cm	10m
FMS25	Green	Medium Limb	5cm*	10 – 25cm	10m
FMS35	Blue	Large Limb	7.5cm*	20 – 45cm	10m
FMS45	Yellow	Extra Large Limb, Head	10.75cm*	35 – 65cm	10m
FMS55	Beige	Adult Trunk	17.5cm*	50 – 120cm	10m

\*Unstretched Width



Fibre content: 92% Viscose, 5% Elastodiene, 3% Polyamide

### Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)

ARTG 222195

Version 2 / Created 18/5/2020