



sutherlandmedical

Introducing Aluminium **Finger** **Splint:**



SPLINTS AND IMMOBILISATION



sutherlandmedical

Aluminium Finger Splints

Sutherland Medical Aluminium Finger Splints can be bent to the desired shape to suit the digit. Prior to application ensure that any open wounds are cleaned and dressed. Select the appropriate sized finger splint and secure into position.

Features:

- 100% Australian Made
- 1.2mm Aluminium strip thickness
- No sharp edges – finger splints are rounded and bevelled
- Padded with a soft, breathable foam
- Non-autoclavable & designed as single patient use only

| Code | Description | Width | Length | Qty |
|------------|-------------------------|-------|--------|-----------|
| FS12480-15 | Aluminium Finger Splint | 12mm | 480mm | 15 / pack |
| FS19480-15 | Aluminium Finger Splint | 19mm | 480mm | 15 / pack |
| FS25480-15 | Aluminium Finger Splint | 25mm | 480mm | 15 / pack |
| FS32480-15 | Aluminium Finger Splint | 32mm | 480mm | 15 / pack |

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au

ARTG 150501

Version 2 / Created 18/5/2020