



sutherlandmedical

# Introducing **Keep Dri** Dressings:



WOUND CARE



sutherlandmedical

# Keep Dri Dressings

**Keep Dri Dressings** are designed to protect plaster of paris and fibre glass casts and other sensitive areas from moisture during showering or bathing.

The **Keep Dri Dressing** rectangles are ideal for application over wound dressings, surgical sites, and stitches.

**Keep Dri Dressings** are 100% Australian Made

Code	Type	Size	Application	Qty
KDB-Y	Bag	60 x 30cm	Lower Leg, Full Arm	25
KDB-B	Bag	100 x 38cm	Full Leg	25
KDR-P	Rectangle	18 x 13cm	Catheters, Key Hole Surgery	25
KDR-Y	Rectangle	23 x 17cm	Small Sized Wound	25
KDR-B	Rectangle	26 x 23cm	Medium Sized Wound, Drainage Sites	25
KDR-G	Rectangle	42 x 23cm	Large Sized Wounds, Hip Replacements	25
KDB-Y5	Bag	100 x 38cm	Lower Leg, Full Arm	5
KDB-B5	Bag	18 x 13cm	Full Leg	5
KDR-Y5	Rectangle	23 x 17cm	Catheters, Key Hole Surgery	5
KDR-B5	Rectangle	26 x 23cm	Small Sized Wound	5
KDR-G5	Rectangle	42 x 23cm	Large Sized Wounds, Hip Replacements	5

## Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)

ARTG 119590

Version 3 / Created 18/5/2020