



sutherlandmedical

Introducing **Relevo™**: The Super Absorbent Dressing.



WOUND CARE

Relevo™

The super absorbent dressing for heavily exuding wounds

Relevo™ is a sterile dressing which absorbs 2-3 times more than traditional absorbent dressings.

The design allows the dressing to be applied either side facing the wound.

Relevo™ absorbs under compression and binds the fluid inside the dressing without leaking.

Features:

- Contributes to fewer changes, time saving and better patient comfort
- Absorbs vertically and can reduce the risk of maceration on healthy skin area
- Contains a combination of cellulose fibres and superabsorbent polymers (SAP) with a polypropylene cover
- Does not contain or made with natural rubber latex.
- Low profile and easy handling

Code	Description	Absorption capacity		Qty
		(Tap water)	(Saline Solution 0.9%)	
RLF50011/10	Relevo™ Dressing 10cm x 10cm	120ml	50ml	10
RLF50011/50	Relevo™ Dressing 10cm x 10cm	120ml	50ml	50
RLF50021	Relevo™ Dressing 10cm x 20cm	300ml	120ml	50
RLF50041	Relevo™ Dressing 20cm x 20cm	600ml	240ml	50



Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au