



sutherlandmedical

Introducing **SMT** Tubular Latex Compression Bandage:



COMPRESSION BANDAGES



sutherlandmedical

SMT Tubular Latex Compression Bandage

SMT is ideally designed for the treatment of orthopaedic, vascular conditions and soft tissue injuries.

SMT is used for:

- Treatment of Venous insufficiency requiring low – medium compression
- Control of Oedema / Swelling of the limbs
- Soft tissue, joint injuries or rib injuries
- Application following the removal of casts

Features:

- **100% Australian Made**
- Latex based thread
- Product treated with Unique Low Fray Formula
- 12 metre rolls
- Easy to use Applicators available

Measuring Guide for Application

Code	Colour	Size	Application	Circumference**
SMT-A	Natural	A	Infant Limbs	7 - 11cm*
SMT-B	Natural	B	Small Limbs	12 - 14cm*
SMT-C	Natural	C	Adult Limbs	15 - 25cm*
SMT-D	Natural	D	Large Limbs	26 - 35cm*
SMT-E	Natural	E	Adult Legs	36 - 44cm*
SMT-F	Natural	F	Large Legs	45 - 50cm*
SMT-G	Natural	G	Large Thighs	51 - 60cm*
SMT-AB	Flesh	A	Infant Limbs	7 - 11cm*
SMT-BB	Flesh	B	Small Limbs	12 - 14cm*
SMT-CB	Flesh	C	Adult Limbs	15 - 25cm*
SMT-DB	Flesh	D	Large Limbs	26 - 35cm*
SMT-EB	Flesh	E	Adult Legs	36 - 44cm*
SMT-FB	Flesh	F	Large Legs	45 - 50cm*
SMT-GB	Flesh	G	Large Thighs	51 - 60cm*

LATEX

*When applying to limbs, measure to fit narrowest part. When applying to other areas of the body, measure the affected area

**Guide is for Medium Compression, for Low Compression select next size up.

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au

ARTG 192869

Version 3 / Created 18/5/2020