



sutherlandmedical

Introducing **Trio Siltac**[®]:
silicone ostomy seal for skin
friendly protection.



Feel the difference, experience the benefit of
our advanced silicone technology.....

OSTOMY CARE



sutherlandmedical

Siltac[®] silicone ostomy seal

Advanced Silicone Technology

To date, you have had a limited choice in ostomy seals. All current products comprise of a hydrocolloid formulation. All current products suffer from the same set of problems.

Trio Siltac[®] is a soft silicone ostomy seal that works completely differently from hydrocolloid products. It does not absorb moisture or body waste for the first time you can have a clean and hygienic seal.

Trio Siltac[®] seals have a natural tendency to return to their original shape when stretched. This means you simply stretch a **Trio Siltac[®]** seal around your stoma, and shape to form a perfect custom fit.

Trio Siltac[®] removes in one piece leaving no residue - your pouch change regime can now be quicker, easier and cleaner

Features:

- Maintains integrity and shape to provide longer lasting seal
- More hygienic than traditional hydrocolloid products - does not absorb faeces or urine
- Allows the skin to breathe
- Reduced risk of leakage - extends wear time of your ostomy product
- Stretch around your stoma and shape to create the perfect custom fit
- **Trio Siltac[®]** conforms to the shape of your stoma
- Promotes comfort due to high level of flexibility
- Instant adhesion - better seal formed
- One piece removal - leaves no residue
- Atraumatic removal - less skin irritation



Code	Description	Qty	SAS Code
TR1020	Siltac [®] Silicone Ostomy Seal, Size 1, 20-28mm	30	80055K
TR1028	Siltac [®] Silicone Ostomy Seal, Size 2, 28-35mm	30	80055K
TR1035	Siltac [®] Silicone Ostomy Seal, Size 1, 35-44mm	30	80055K

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au

ARTG No 278399

Version 2 / Created 18/5/2020