



sutherlandmedical

Introducing **Tubular - Band™**:
A clever alternative to
conventional dressing retention.



RETENTION BANDAGES

Tubular - Band™

Tubular Retention Bandage

Tubular - Band™ is ideally designed as a low cost, easy to use alternative to crepe bandage and adhesive applications

Tubular - Band™ is used for:

- Wound dressing retention
- Isolation of I.V. Lines and Tubing
- Application to any area where Skin Protection is required

Features:

- 100% Australian Made
- Extremely lightweight bandage
- Provides effective retention for every part of the body
- Quick and easy to apply
- Fully washable
- Extremely cost competitive
- Available in 5 easy to identify colour coded sizes

Measuring Guide for Application

Code	Colour	Limb Circ	Application	Flat Measure	Roll
SM900	Red	9 – 20cm	Small Limbs	3.5cm*	10m
SM901	Green	14 – 25cm	Medium Limbs	5cm*	10m
SM902	Blue	24 – 36cm	Large Limbs	7.5cm*	10m
SM903	Yellow	35 – 65cm	Head or Small Trunk	10.75cm*	10m
SM904	Beige	64 – 100cm	Large Trunk	17.5cm*	10m



*Size= unstretched width.

Fibre content: Viscose, Nylon, Elastodiene

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au