




sutherlandmedical

Introducing Tubular - Rib™: A natural alternative for skin protection.



Insert for all rolls



PROTECTION BANDAGES

Tubular - Rib™

Tubular Retention Bandage

Tubular - Rib™ is ideally designed to provide a protective layer in direct contact with the skin.

Tubular - Rib™ is used for:

- Protection of the skin prior to the application of fibre-glass and plaster of paris casts.
- Protection of the skin under all types of compression bandaging.

Features:

- 100% Australian Made
- Cotton 100%

Measuring Guide for Application

Code	Flat Measure	Application	Limb Circ	Roll
SM500	5cm*	Small Limb	15 – 22cm	20m
SM501	7.5cm*	Medium Limb	20 – 32cm	20m
SM502	10cm*	Large Limb	30 – 47cm	20m
SM503	15cm*	Extra Large Limb	45 – 62cm	20m
SM504	20cm*	Small Trunk	60 – 78cm	20m
SM505	25cm*	Large Trunk	75 – 90cm	20m



*Size= unstretched width.

Fibre content: Cotton Elastodiene, Nylon

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au