



sutherlandmedical

# Introducing **Nitto** Yu-ki Ban™ & Xtrata™ Perme Roll: Surgical tapes & transparent film dressings with STRATAGEL™ technology.



WOUND CARE




## Nitto Yu-ki Ban™ & Xtrata™ Perme Roll

### Stratagel™ adhesive technology, protecting the skin from both irritation and pain on removal

Above all, Gentle to the skin  
Stratum Corneum Protection System

**Conventional Product**  
Skin damage caused by removal of stratum corneum

**XTRATA**  
Very little Stratum corneum is removed



99% of stratum corneum is removed  
Conventional Product A

40% of stratum corneum is removed  
Xtrata Perme-roll

#### Features:

##### XTRATA PERME-ROLL

- Acts as a barrier to external bacteria
- Superior moisture permeability
- Excellent adhesion with low removal of stratum corneum
- Superb flexibility and elasticity to joints and other flexible areas
- Waterproof
- Transparency allows observation of the application site without removal
- Latex Free

##### YU-KI BAN

- Low irritation tape
- Allows repeated application
- Latex Free



#### Code

#### Description

#### Qty

RL3251	Yu-ki Ban™ Non-woven (water resistant) cloth 12mm x 10m	24 rolls
RL3252	Yu-ki Ban™ Non-woven (water resistant) cloth 25mm x 10m	12 rolls
RL3255	Yu-ki Ban™ Non-woven (water resistant) cloth 50mm x 10m	6 rolls
RL3521	Yu-ki Ban™ Alpha (waterproof) polyolefin film 12mm x 7m	24 rolls
RL3522	Yu-ki Ban™ Alpha (waterproof) polyolefin film 25mm x 7m	12 rolls
RL3525	Yu-ki Ban™ Alpha (waterproof) polyolefin film 50mm x 7m	6 rolls
RLH24R05	Xtrata™ Perme Roll (transparent film dressing) 5cm x 10m	1
RLH24R10	Xtrata™ Perme Roll (transparent film dressing) 10cm x 10m	1

#### Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)

Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)