



sutherlandmedical

Introducing CollarForm™ Arm Sling:



SPLINTS AND IMMOBILISATION

CollarForm™

CollarForm™ is a strong, stretch resistant sling made of a 5cm wide strip of dense foam covered with a pure cotton stockinette.

CollarForm™ is ideally designed to immobilise the clavicle.

CollarForm™ is 100% Australian Made

Code	Width	Length	Packaging
SM000	5cm	10m	Economy Roll
SM001	5cm	2 x 6m	Dispenser Pack



Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189
Ph: 1300 664 027
Email: orders@sutherlandmedical.com.au
Web Site: www.sutherlandmedical.com.au