



sutherlandmedical

# Introducing the **Focus™** range of stethoscopes:



Figure 1 Classic Black



Figure 2 Classic Red



Figure 3 Classic Green



Figure 4 Deluxe



Figure 5 Cardio

DIAGNOSTICS



sutherlandmedical

# FOCUS™ Stethoscopes

The **FOCUS™** range of stethoscopes is a professional series of Stethoscopes designed to provide ideal aural quality within the hospital or specialist environment.

## Features:

- Internal Spring Binaurals
- Matching tubing
- Non chill ring

The **FOCUS™** Classic stethoscope features a new design chrome plated zinc alloy chest piece and a spare pair of ear pieces.

The **FOCUS™** Deluxe stethoscope features a premium quality stainless steel chest piece providing superior quality sound transmission with a high quality diaphragm and a spare pair of ear pieces

The **FOCUS™** Cardiology stethoscope features a premium quality stainless steel chest piece providing superior quality sound transmission with a high quality diaphragm and is supplied with a spare diaphragm and a spare pair of ear pieces

The **FOCUS™** Paediatric stethoscope features a premium quality stainless steel chest piece providing superior quality sound transmission and a spare pair of ear pieces  
(Nylon 2%)

Code	Product	Colour	Carton
SMS0010	Classic	Black	50
SMS0020	Classic	Red	50
SMS0030	Classic	Green	50
SMS0040	Deluxe	Black	50
SMS0050	Cardiology	Black	50
SMS0060	Paediatric	Black	50
SMS0090	Ear Pieces	White	Pair



## Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189  
Ph: 1300 664 027  
Email: [orders@sutherlandmedical.com.au](mailto:orders@sutherlandmedical.com.au)  
Web Site: [www.sutherlandmedical.com.au](http://www.sutherlandmedical.com.au)