



sutherlandmedical

Introducing **Tubular - Gauze™**:
A low cost alternative for skin
protection.



PROTECTION BANDAGES



sutherlandmedical

Tubular - Gauze™

Tubular Gauze Bandage

Tubular - Gauze™ is ideally designed to provide a protective layer in direct contact with the skin.

Tubular - Gauze™ is used for:

- Protecting clothing and bedding from soiling by ointments and creams
- Direct wound dressing for fingers and toes
- Providing a bulk protective dressing for fingers
- Providing a protective layer prior to application of adhesive strapping

Features:

- 100% Australian Made
- Cotton 100%
- Extremely light weight bandage
- Ideal for podiatry application
- Applicators are available for fingers and toes

Measuring Guide for Application

Code	Size	Limb Circ	Application	Flat Measure	Roll
SM600	00*	2.5 – 5cm	Small Fingers & Toes	0.8 cm	20m
SM601	01*	5.5 – 7.5cm	Fingers & Toes	1.8 cm	20m
SM602	12*	7.5 – 12cm	Bulk Finger Dressings	2.5 cm	20m
SM603	34*	12.5 – 20cm	Child Limb	4 cm	20m
SM604	56*	20.5 – 25.5cm	Adult Limbs	6.5 cm	20m
SM605	78*	25 – 37cm	Large Limbs	8.5 cm	20m

*Flat Measure = unstretched width.



Fibre content: Cotton

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au

ARTG 133290

Version 4 / Created 16/10/2020