



sutherlandmedical

Introducing **Tubular - Net™**:
A low cost alternative to
conventional dressing retention.



RETENTION BANDAGES

Tubular - Net™

Tubular Net Retention Bandage

Tubular - Net™ is ideally designed as a low cost, easy to use open weave retention bandage.

Tubular - Net™ is used for:

- Wound dressing retention
- Isolation of I.V. Lines and Tubing
- Any retention application where minimum skin contact is required

Features:

- 100% Australian Made
- Extremely lightweight bandage
- Provides effective retention for every part of the body
- Ideal for Burns Units and Neonatal Intensive Care Units
- Constructed Free of Knots eliminating risk of pressure damage
- Applicators are available for fingers and toes

Measuring Guide for Application

Code	Size	Limb Circ	Application	Flat Measure	Roll
SM300	0	2 – 6cm	Fingers & Toes	0.8cm*	25m
SM301	1	6 – 9cm	Lower Arm	1.9cm*	25m
SM302	2	9 – 11cm	Hand or Foot	2.7cm*	25m
SM303	3	11 -13cm	Arm or Leg	3.5cm*	25m
SM304	4	13 – 15cm	Knee or head	4.0cm*	25m
SM305	5	14 – 26cm	Shoulders or Head	4.4cm*	25m
SM306	6	26 – 33cm	Adult Chest (blue line)	7.8cm*	25m
SM307	7	30 – 38cm	Large Chest (yellow line)	9.6cm*	25m
SM308	8	32 – 45cm	Extra Large Chest (red line)	10.4cm*	25m
SMAPP			Finger applicator		



*Flat measure= unstretched width.
All rolls are stretched length (1m unstretched = 2.5m stretched)

Fibre content: Nylon, Elastic

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au

Web Site: www.sutherlandmedical.com.au

ARTG 133287

Version 3 / Created 1/6/2021