


KIS™ AMORPHOUS HYDROGEL

Gentle amorphous care for facilitating the wound healing process.



KIS Amorphous Hydrogel is a colourless transparent gel. The gel can gently rehydrate necrotic tissue and facilitate autolytic debridement. It contains water and propylene glycol which allows it to donate moisture to dry wounds. The gel creates a moist wound healing environment.

- Carboxymethylcellulose and calcium chloride.
- Amorphous Hydrogel contains water and propylene glycol.
- Gently rehydrates necrotic tissue.
- Easy to use, single use only.
- Facilitates autolytic debridement.
- Creates a moist wound healing environment.



HYDRATION

Indications:

KIS Amorphous Hydrogel is used to rehydrate necrotic tissue and facilitate autolytic debridement, while not damaging fragile granulation tissue:

- Leg ulcers
- Pressure Injuries
- Diabetic foot ulcers
- Wound Cavities

Contraindications:

- For external use only.
- Cannot be used to heal burns.
- Do not use if a known sensitivity to any of the components.

Benefits:


- Cooling and soothing
- Easy to apply




Moisture balance

Codes	Sizes
SMAH015B	15g Tube
SMAH025	25g Tube


APPLICATION INSTRUCTIONS:



1 Snap off the cap from tube.



2 Apply Hydrogel to cover wound bed.



3 Secure with secondary dressing if needed.